GYMSTICK™

Cair0


GYMSTICK™
MASSAGE GUN
CAIRO PRO

Cair0

USER MANUAL
Percussion Technology

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MASSAGE GUN CAIRO PRO

DISCLAIMER

The manufacturer Gymstick International Oy have no responsibility to any persons or entity for any loss, damage sequence alleges to have happened directly or indirectly as a consequent of this product.

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Gymstick International Oy warrants this product to be free of defects in material or workmanship for a period of one (1) year from the date of the original retail purchase.

This warranty is available only to the original purchaser, and only if purchased from Gymstick distributor or reseller. If your product is found to be defective during the warranty period, Gymstick will, at its sole option, either repair or replace the product with a new or factory reconditioned or substantially equivalent model. Gymstick will not replace or refund cover damage of the body due to consumer's negligence and carelessness.

The guarantee does not cover: Damage effected by outer force, Intervention by unauthorized parties, Incorrect handling of the product and Non-compliance of the operating instructions.

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- Operating time: about 5 hours
- Charging time: about 2 hours
- It may change depending on the usage environment

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

Mitat	18,5 x 27 x 7,6 cm
Paino	967 g
Nimellisteho	35W harjaton moottori (maks. 80W)
Akku	24V / 2600mAh (Li-ion)

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SPECIFICATIONS & FUNCTIONS

OPERATION AND FUNCTION

- There is ON and OFF power switch at the bottom of the battery handle
- The charging indicator will show the remaining power percentage when power ON.
- The device will shut down automatically when continuously using 30 min. (to protect the battery life).
- When not in use for more than 3 min. the device will turn to sleep mode automatically.
- Touch screen displays the cumulative number of beat times (every 100 beatings show 1 digit).

NOTE: Switch to "OFF" position when not using it and place the device properly.

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
Press "+" to increase intensity / to start massage
Press "-" to decrease intensity / till the device stop working

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CHARGING INSTRUCTION

- The battery is fully charged before leaving manufacturer.
- The power indicator light will flash in RED when the electricity remains 10%. When the remaining power is 0%, the power indicator light starts to flash and the machine will automatically shut down.
- It will take about 3 hours to fully charge the battery.

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
ASSEMBLE INSTRUCTION

1. Insert the battery handle to the main body. The battery is fully connected when you hear "click" sound.

Dissembling: Press the button on main body to detach the battery handle from the main body.
(The battery may feel tight and difficult to remove. This is not a product defect but to minimize the gap between the battery and the main body)

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

3. Replace the massage head in time according to different massage parts.

(if installation method is wrong, the massage head may fall off during operation. Please make sure the massage head is tight before using.)

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

! When install the **U-head** massage head. Please make sure the U-head is pushed directly into the interface groove of the massage device and align the mortise and the tenon joint. Do not spin the U-head to avoid abrasion of the rubber band.

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FUNCTION OF MASSAGE HEADS

Big ball head:

Suitable for whole body, especially for large muscle groups such as waist, hip and leg.

U-head:

Suitable for shoulder and spine part or enhance for key areas.

Air cushion:

Suitable for whole body, especially for biceps muscles and leg.

Flat head:

Suitable for whole body, especially for waist, chest, inner and outer side of leg.

Bullet head:

Suitable for feet, hands, Achilles' Tendon and back muscles.

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EXERCISE MANUAL

PLEASE READ THE FOLLOWING CAREFULLY BEFORE ASSEMBLING AND OPERATING THIS MACHINE

1. This device cannot be used as a self-diagnosis treatment. If it needs diagnosis or treatment, it should be handled by professional medical personnel.
2. If the device is used or close to children, physical and mental dysfunction and the disabled, it is necessary to have guidance from others and be responsible for its safety.
3. Do not operate massage device in a humid or explosive atmospheres, such as in the presence of flammable liquids, gases or dust. Massage devices create sparks which may ignite the dust or fumes.
4. This machine should not be exposed to the rain or immersed in water. The main body, vent, battery handle and adaptor, please keep away from all kinds of liquid.
5. Do not crush the product and battery and do not place the device near the electrical heater, near the stove, humid atmosphere.
6. Do not use this device in humid or high temperature environment. Do not use the device when the room temperature is higher than 40 degrees.
7. Please stop using the device when the machines falls down to the ground, into the water and have defects.
8. Do not drop any small items from the device cooling window or opening into the device interior.
9. Do not place or use the device under a pillow or blanket.
10. It may cause discomfort, when the user overuses the device. The device will shut down automatically when continuously using 30 minutes.
11. Forbidden to use on the head, abdominal cavity, joints or bones.
12. Do not overuse the device to avoid hurting.
13. Dress properly. Do not wear loose clothing or jewelry. Keep hair, clothing and gloves away from moving parts.
14. Do not place the device under direct sunlight.
15. Do not continuously use more than an hour. Please rest for 30 minutes then continue to use it again.

If you have any of the following conditions, please consult your doctor before using!

Hypertension, heart disease, cerebral arteriosclerosis, malignant tumor, fractures or osteoporosis, pacemaker or implanted medical device, you are pregnant or you suffer from any acute disease.

Please stop using the machine when you feel uncomfortable.

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- If the device is used or close to children, physical and mental dysfunction and the disabled, it is necessary to have guidance from others and be responsible for its safety.
- Recharge only with the charger specified by the manufacturer.
- Do not operate the control switch with wet hands.
- Always turn off the device when stop using it.
- Do not fall asleep while the device is in operation.
- Do not use a massage device while you are tired or under the influence of drugs, alcohol or medication.

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How to: Swing the neck naturally. Massage the shoulder and neck from top to bottom and back and forth.

Time: 20 to 30 seconds for each part. The overall massage time is recommended to be 2 to 5 minutes.

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
INTENSITY: 1400 - 2400 rpm

Deltoid muscle / Hartialihäs / Deltamuskeln

How to: Lift your arm up and stretch it, swing it naturally from top to bottom. Place the massage device on the front of the deltoid muscle for massage.

Time: 20 to 30 seconds for each part. The overall massage time is recommended to be 2 to 5 minutes.

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INTENSITY: 1400 - 2400 rpm

Biceps / Haislihas / Biceps

How to: Straighten your arm forward. Make a fist. Fold your forearm up and massage along the inside of the upper arm.

Time: 20 to 30 seconds for each part. The overall massage time is recommended to be 2 to 5 minutes.

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INTENSITY: 1400 - 1800 rpm

Rhomboid muscles / Suunnikaslihas / Rhomboideus muskeln

How to: Fold or stretch your arms back and forth, massage the muscles between the shoulder and the spine.

Time: 20 to 30 seconds for each part. The overall massage time is recommended to be 2 to 5 minutes.

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How to: Bend the body forward slightly, hold the massage device backhanded and massage the waist.

Time: 20 to 30 seconds for each part. The overall massage time is recommended to be 2 to 5 minutes.

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How to: Slide the device naturally to massage from hips to leg by clockwise from 2 o'clock to 9 o'clock. Move the device back and forth.

Time: 20 to 30 seconds for each part. The overall massage time is recommended to be 2 to 5 minutes.

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How to: Massage on the inner side of the thigh. Massage back and forth.
Time: 20 to 30 seconds for each part. The overall massage time is recommended to be 2 to 5 minutes.

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INTENSITY: 1400 - 2800 rpm

Gastrocnemius / Pohjelihak / Vadmuskeln


How to: Slide the massage device along the back of the calf toward the ankle for massage.
Time: 20 to 30 seconds for each part. The overall massage time is recommended to be 2 to 5 minutes.

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INTENSITY: 1400 - 2400 rpm

Plantar muscles / Kantalihakset / Plantarmusklerna

How to: Place the massage device on the soles of the feet and massage with the forward and backward movements of the toes. Stretch your toes forward and back.

Time: 20 to 30 seconds for each part. The overall massage time is recommended to be 2 to 5 minutes.

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INTENSITY: 1400 - 1800 rpm

Pectoralis / Rintalihas / Pectoralis muskeln

How to: Stretch the arm to the position of 9 and 11 o'clock and massage the pectoralis muscles gently.

Time: 20 to 30 seconds for each part. The overall massage time is recommended to be 2 to 5 minutes.

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INTENSITY: 1400 - 1800 rpm

Time: 20 to 30 seconds for each part. The overall massage time 2 to 5 minutes.

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Troubleshooting

If you have problems using this device, please refer to the following tips and solutions. If the device is still in a fault condition, please contact with our customer service department (info@gymstick.com).

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The device suddenly stops during operation.	
Battery temperature too high.	Detach the battery handle from the main body and cool it for about 30 minutes.
The device has been in use for over 30 min.	Switch to "OFF" and then switch to "ON" again. Then restart using the device.
The device is not working.	
The battery is running out of electricity.	Please charge the device. Please note that the device won't operate during the charging.
The device is not turned on.	Please switch to "ON".
The device is not reacting even though the switch is turned "ON" again.	
Depending on the use environment, it is difficult to determine the cause.	Charge the battery again and check if the indicator of the battery is "ON". After confirming the battery level, unplug the adaptor and operate the machine.
The percentage of the battery consumption show on the bottom, but without reaction when touch the screen.	
It may be that the battery is not connected properly to the main body.	Detach the battery handle first and fully connect it till you hear "CLICK" sound.

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Maintenance:

- Please use a dry cloth to clean the surface of the device. Never use water or chemical cleaners.
- Check for a breakage and any other condition that may affect the massage device's operation.
- When not in use, place the device in the travel case and keep it in a safe, dry, ventilated, cool place. Do not place the device in direct sunlight or high temperatures.

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Do not dispose of the device in the normal household waste, but in accordance with the Directive and local regulations. More information on this is available from your municipal waste collector.

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