
W

Svenska Hälsobänken
The Swedish Back Care System
Die Schwedische Gesundheitsbank
La table thérapeutique suédoise
El banco de salud sueco

Bruksanvisning
Operating instructions
Anwendungsanweisungen
Mode d'emploi
Modo de empleo

Model CN-A, CN-B

Made in China for Mastercare AB
PAT. / DES. PAT, WORLD WIDE

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Congratulations on your decision!
Welcome to a healthier life with Mastercare - The Swedish Back Care System

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This instruction manual should always be kept available next to
Mastercare - The Swedish Back Care System.
Please read the manual carefully and follow the instructions.

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Pictures and descriptions in this publication were correct at time of printing. We reserve the right to make changes due to a continuous development and improvement of our products.

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Mastercare - Svenska Hälsobänken

Instruction Manual

MASTERCARE – THE SWEDISH BACK CARE SYSTEM

Mastercare – The Swedish Back Care System is a tool, which helps you to prevent and relieve pain from back, neck, shoulders and knees, as well as wearing out of joints.

The construction of Mastercare – The Swedish Back Care System with a movable back rest, transforms your body weight into traction power, which relieves discs and joints. Simultaneously, movements on the table involves strain of muscles, muscular attachments, tendons, ligaments and cruciate ligaments. At the beginning this will cause slight stiffness and ache, a sign that the parts are too short or too weak. These parts need to be built up and therefore it is important not to stop training, even if/when you will be stiff after training. Daily use of Mastercare – The Swedish Back Care System can give you a healthier life.

The table is designed for use by one person at a time. When you start your training for the first time, use the table in the presence of a person, who is familiar with the functions of the table.

Be sensitive to the signals from your body. Consult with your doctor or physiotherapist, if you start aching, feel uncomfortable or notice any unusual changes that you think might prevent you from using the table. Don't exercise on the table immediately after a meal or if you have a temperature, a cold or severe head- or earache.

If you are uncertain about the state of your health or suffer from any of the following conditions, **ALWAYS** consult your doctor before starting a regular use of the table:

Glaucoma, detached retina, conjunctivitis, uncontrolled hypertension (high blood pressure), impaired pulmonary functions, arrhythmia, carotid artery stenosis, severe vascular disease, aneurysm, tendency towards fainting or dizzy spells, cerebral sclerosis, recent stroke or transient ischemic attack, history of brain lesions, other heart and circulatory disorders for which you are being treated, severe back problems or spinal injuries, bone weakness (osteoporosis), recent unhealed fractures, medullary pins, artificial hip joints or other orthopedic supports surgically implanted, middle ear infection, recently had major surgery, using special medication, ventral or hiatal hernia, tendency towards bleeding, extreme obesity, invertaphobia (fear of inversion itself).

Using the table

Before using the table on your own, you should acquire the proper know-how.

Check that:

- **the table is mounted with stretched cords (see assembly instructions 1B)**
- **it stands on a level, horizontal surface**
- **all nuts and bolts are securely fastened**
- **the foot lock is working perfectly**

Mastercare – The Swedish Back Care System is designed for persons with a height of 135-215 cm and a maximum weight of 150 kg.

**Note! Always be careful so you don't get your fingers jammed.
Always keep your hands on the handles when the table is in tilting motion.
Don't leave children using the table without supervision.**

Find out how the safety catches work and make sure that at least one of the catches **always** is in function. The other catch can be disengaged (picture 1), if you want to operate the table with just one hand. If both catches are disengaged, the table will swing without safety stop and you may be injured.

In order to operate the table with the best balance, you should adjust the foot rest to the length corresponding to the length of your body (picture 2). If this position doesn't give good balance, adjust the foot rest to the "length" that will give you the best balance.

Wiggle slightly with the table when you loosen the safety catches. If it works stiffly, the length is not properly adjusted. When lying in prone position (on your stomach) on the table, you should always be assisted by a person who is familiar with the functions of the table.

Note! When lying in prone position, and are returning to the initial position, do this carefully or you might injure your toes. By length adjustment between 185-215 cm (model A), toe protection should be used (extra accessory). Model B for supine position only.

Before stepping on to the foot rest, always check that the catch for the length adjustment is in locked position (picture 2).

When you are about to use Mastercare – The Swedish Back Care System for the first time, it is appropriate to make a "handling exercise" in order to feel confident with the functions of the table. Here is how:

- Take off your shoes.
- Step on to the foot rest with your back against the table and your heels against the foot rolls (picture 3), or facing the table, when you use it as a massage table or work in prone position.
- Tighten the foot rolls to your insteps.
- Lean backwards against the back rest.
- Press down the locking devices on both sides of the stand, and press yourself slowly backwards (picture 4).
- When the back rest starts tilting backwards, let go of both catches, but keep your hands pressed against the handles, until the catches will lock in the next position.
- Continue to horizontal position (picture 5), where you should make a halt, both at the beginning and the end of the exercises.
- Press down the catches and continue to 15° inclination (picture 6), which is the first position with your head lower than your feet.
- When the back rest is locked, place your hands on your stomach and bend/stretch both legs, so that you can feel the back rest move on its rolls.
- Press down both catches and return, **with a break in the horizontal position where you rest**, to the initial position (picture 19).

Repeat this exercise until you feel comfortable using Mastercare – The Swedish Back Care System. Now this you are ready to start using the table.

A training programme

Now bear in mind that it is your body that you are working with. Follow the instructions you have received, work quietly and carefully, and be sensitive to the signals of your body. For best result it is necessary that you can relax and take some deep breaths between the different parts of your movements/exercises.

“Repeat” will return in many parts of the instruction. It is for you to decide if you want to repeat the exercises one time or more. An easy start is recommended. Increase your muscular tension when you repeat an exercise!

Think positively. Start now!

The numbers to the left are referring to the pictures on page 3.

2. Adjust the length for correct balance and lock the catch.
3. Step on to the foot rest and tighten the foot rolls to your insteps.
4. Press down the locking devices and press yourself backwards to the horizontal position.
5. Take some deep breaths. Breathe in, hold your breath, breathe out slowly and try to relax. Stretch out with your hands pressed against the handles. Repeat!
6. Continue to 15° inclination, which is the first position with your head lower than your feet. The following exercises are to be performed in this position.
7. Stretch out from top to toe, with your arms down over your head. Then fold your hands on your stomach. Take some deep breaths. Relax.
If you feel comfortable after some deep breaths, you are ready to start with the movements. If you don't feel comfortable, return to horizontal position. Relax. After a short break, try again. If you still don't feel comfortable, finish the exercise (according to points 19 - 22), and step down carefully from the table.
8. Press the upper side of your feet against the foot rolls with straight knees. Hold for a moment. Relax. Take some deep breaths. Repeat.
9. Lift one knee as high as you can. Keep the other leg straight. Hold for a moment. Slowly bring your leg down. Lift the other knee. Relax. Take some deep breaths. Repeat.
10. Lift your seat from the back rest and pinch. Hold for a moment. Lie down. Relax. Take some deep breaths. Repeat.
11. Cross your arms over your chest. Roll up on one shoulder and stretch your back. Hold for a moment. Lie down for a moment. Roll up on the other shoulder. Hold for a moment, then down on your back. Relax. Take some deep breaths. Repeat.
12. Stretch your arms and hold the palms of your hands slightly pressed together over your stomach. Slowly bend your back upwards. (*Note! No sit-up!*) Hold for a moment. Lie down. Carry out the movement with your hands stretched at first over the right handle, then over the left one. Relax. Take some deep breaths. Repeat.
13. Press your head against the back rest and lift your shoulder blades. Hold for a moment. Lie down. Relax. Take some deep breaths. Repeat.
14. Stretch out your arms. Twist your arms as far as you can towards both sides. Spread your fingers and clench your hands alternately while twisting. Move your outstretched arms between five different positions and carry out the movements. Relax. Take some deep breaths. Repeat.
15. Put your hands on your shoulders. Roll your shoulders slowly forwards and backwards. Relax. Take some deep breaths. Repeat.
16. Clasp one wrist and pull the arm over your head in the length direction of the table. Do the same with the other arm. Relax. Take some deep breaths. Repeat.
17. Put your hands on your stomach. Turn your head as far as possible towards each side. Relax. Take some deep breaths. Repeat.

18. Fold your hands under the back of your neck and press your chin against the chest, without lifting your shoulders. Hold for a moment.
Relax. Take some deep breaths. Repeat.

Now you have **one** exercise programme.

Start with caution. In the beginning, "hold for a moment", means a few seconds. After you have been training for some weeks, it may include 20 - 30 seconds. Feel how your body reacts. Intensify the exercises gradually.

Feel free to talk to your back specialist about an exercise programme just for you.

After the exercise it is important, **very important**, to resume the initial position, without anything in your back getting jammed.

Do like this:

19. Pull the back rest up to horizontal position.

Relax. Take some deep breaths. Rest a while (picture 5).

20. Pull the back rest to the next position. Bend one leg at a time, a few times, to increase the blood circulation. Stand/lie straight with your arms along the sides. Pull one shoulder upwards, the other one down, alternately a few times.

21. Pull the back rest to the initial position. Bend your knees. Repeat the previous movement (picture 20).

22. Release the foot lock and cautiously step out of the table.

After some time, when you are familiar both with Mastercare – The Swedish Back Care System and the reactions of your body, you may try some of the exercises in 30° inclination. For a perfect result, total relaxation is necessary, and it is often most easily obtained in 15° inclination.

We congratulate you on your decision to help your body to a healthier life!

Mastercare – The Swedish Back Care System

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