
PT Board

Manual
User handbook

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CONGRATULATIONS on purchasing PT BOARD from Sableline Healthcare! To get the most out of this product, please read these instructions carefully before using.

We hope you will enjoy using your (P)ersonal (T)rainer BOARD.

SAFETY INSTRUCTIONS

Please retain these user instructions.

- Your PT BOARD complies with the recognised principles of technology and current safety regulations.
- All necessary repairs must be carried out by qualified persons. Unauthorised repairs and faulty use may render the warranty invalid. Always contact your retailer, see warranty certificate.
- Never touch the electrical plug with wet hands. If the equipment is going to be out of use for an extended period, make a point of unplugging it from the mains socket. This is partly to save power and also to avoid short-circuiting. Never use damaged plugs or cord. A damaged plug should be replaced by your retailer, see warranty certificate. In the event of malfunction indication, disconnect the unit from the mains socket immediately and contact your retailer. Sableline Sweden AB does not accept any liability for damage resulting from incorrect handling.
- Do not allow the equipment to come into contact with intense heat, direct sunlight or water.
- If the equipment has been stored in a cold place for an extended period, allow it to stand at room temperature for at least one hour before starting up.

Even the most physically fit person should start slowly, as this alternative training system represents an entirely new method of exercising the muscles. The intensity of training can be increased gradually, with the frequency being stepped up as differing exercises are introduced.

The overwhelming opinion is that younger children should not use PT BOARD for exercising.

In order to achieve the best results and avoid injuries, always make sure that you do not have any sharp objects on you or close to where you are exercising. Never place your chest or head directly against PT BOARD. Always follow the enclosed exercise instructions and do not use this alternative training system for treating or alleviating illnesses without first consulting your doctor.

WHEN NOT TO USE PT BOARD

- acute thrombosis, or risk of thrombosis
- pregnancy
- severe brittleness of the bones with fractures
- a recent operation, recent wounds or bone fractures
- acute disc prolapse or inguinal hernia
- untreated high blood pressure
- severe diabetes
- implant, ceramic or metal, less than 6 months old
- pacemaker or other medical electrical implant

REMOTE CONTROL

ON/OFF	For turning PT BOARD ON or OFF
Mode	When PT Board is ON it is possible to choose between 4 different programs: P1, P2, P3, P4 (see automatic programs)
Time	Only in manual mode (1-15 minutes)
H+	Increases speed in amplitude mode (this button does not work in auto mode)
H-	Reduces speed in amplitude mode (this button does not work in auto mode)
L+	Increases speed of lateral movement (this button does not work in auto mode)
L-	Reduces speed of lateral movement (this button does not work in auto mode)

H = amplitude movement, i.e. wave motion

L = lateral movement

H and L can be co-ordinated to produce amplitude movement and lateral movement simultaneously

PT BOARD DISPLAY

ON/OFF	Starts standby mode on PT Board
Mode	Program selection: P1, P2, P3, P4
Speed	H mode (Amplitude) L mode (Lateral movement)
Time	Indicates length of time that PT Board is running
Calories	Indicates average calorie consumption only.

BLUETOOTH - USB

1. Turn on the power and turn on the PT Board. Bluetooth / USB will automatically go into sleep mode.
 2. Open mobile bluetooth, it will automatically connect your mobile Bluetooth device.
 3. Wireless Bluetooth function works up to a distance of 8 meters.
 4. You can not use Bluetooth and USB at the same time.
 5. Button 1: Short press for the previous song, to decrease the volume, hold the button for a while.
 6. Button 2: Short press for the next song, to increase the volume, hold the button for a while.
 7. Button 3: Short press for the next song, to increase the volume, hold the button for a while.
- Support MP1, MP2, MP3

ALTERNATIVE TRAINING WITH PT BOARD

**To achieve the best results, please follow the instructions below.
Do not exceed the recommended times and frequencies.**

physically unfit max 3-5 min 1-4 x exercise sessions/week	average fitness max 10-20 min 1-5 x exercise sessions/week	physically fit max 30 min 1-6 x exercise sessions/week
--	---	---

Dynamic exercise on PT BOARD may also be compared with walking, jogging or running exercises. The position adopted on PT BOARD will determine which amplitude movement is provided and, combined with the speed (1-20), how fast the movement will be.

Example: Position 1, speed 7 = brisk walk.

Positioner

1. walk
2. jogging
3. running

EXERCISE SUMMARY

The recommended exercising on PT Board is extremely dependent on the age and physical condition of the individual user. Since exercising on PT BOARD is highly effective, just brief exercise sessions will be sufficient. A complete exercise routine is provided with each PT BOARD, along with a video illustrating the various exercises involved. Dynamic exercise is very easy on all joints such as hips, knees and ankles. If you have problems with arthritis, PT BOARD is ideal for you, as it will provide you with all the movement you need in a very short time. This will affect and stimulate cartilage formation and increase the production of synovial fluid. Dynamic exercising with PT BOARD involves alternate movements which are adjusted to the body's normal centre of gravity. These trigger a stretch reflex in the muscles which can produce up to 1,000 contractions per second, which means that small micro-muscles throughout the body are also effectively exercised.

WHAT GOALS CAN I ACHIEVE WITH PT BOARD?

- Increased circulation
- Better balance and co-ordination
- Personal training to suit individual ability
- Stronger muscles
- Exercising of micro-muscles
- More effective exercising
- Improved posture
- Stimulation of tissues
- Increased calorie burning

MEDICAL INDICATIONS

- **Defective posture resulting from injury:** Build up your strength and improve your posture.
- **Back pains, various causes:** Improves pain alleviation through muscle relaxation and stretching exercises and co-ordination.
- **Muscle strain, various causes:** Improves muscle performance.
- **Scoliosis:** Improves stability of spine.
- **Excessive muscle tension:** Reduces tension in the muscles.
- **Disc prolapse, non-acute:** Exercises local muscular system.
- **Osteoporosis if not acute and no risk of thrombosis:** Increases bone strength.
- **Joints:** Strengthens joints and improves mobility.
- **Arthritis:** Increases cartilage formation and promotes production of synovial fluid.
- **Parkinson's disease:** Improves co-ordination and posture. (More information can be found under Parkinson).
- **Multiple sclerosis:** Improves posture and increases muscle function, controls muscle tension. (More information can be found under multiple sclerosis).
- **Drop foot syndrome:** Example, disc prolapse (increases muscle performance and improves motor control).

- **Stroke:** Faster mobilisation and rehabilitation. Reduces spasticity, improves power management. (More information can be found under stroke).
- **Balance disorders:** Improves balance and prevents falls.
- **Brachial paralysis:** Example: following an accident, paralysis of nerves in the arm, increases muscle power and performance, improves motor control.
- **Paraplegia:** (Control of tension) Exercising arm and support muscles, relaxation, control of movement, co-ordination, speed, also support for treadmill workout, and incomplete paraplegia. Reduces differences between right and left side when walking, and improves performance.
- **Arterial and venous circulatory disorders (PAOD):** Increases circulation, controls tension in vascular wall muscles. (More information can be found under PAOD).

PRE-PROGRAMMED EXERCISE PROGRAM

Program	Mode	Speed									
		10 minute	9 minute	8 minute	7 minute	6 minute	5 minute	4 minute	3 minute	2 minute	1 minute
P1	H	2	4	6	8	10	12	14	16	18	20
	L	5	8	10	13	15	20	15	13	10	5
P2	H										
	L	2	4	6	8	10	12	14	16	18	20
P3	H	12	10	8	12	10	8	12	10	8	12
	L	8	5	10	5	10	15	5	10	15	10
P4	H	20	18	16	14	12	10	12	10	12	10
	L	8	10	8	10	8	10	15	20	15	20

TECHNICAL DATA

Amplitude	0-18
Frequency	5.6-14 Hz
Speed	1-20
Voltage	AC 220-240V 50/60Hz
Maximum load	120 kg
Power consumption	Max 200 Watt
Weight	23 kg
Dimensions	79 x 46.5 x 14 cm

El-Kretsen →

WARRANTY CONDITIONS

WARRANTY

Provided that PT BOARD is not subjected to direct damage and that the care instructions are followed, a 2- year warranty is provided on both motor and material.

Date of purchase: _____

For easier and faster handling of warranty claims, please use the e-mail address: info@sableline.se.

Describe the fault and give your name, address and telephone number, and attach a copy of the warranty certificate or purchase receipt.

Sableline Sweden AB

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523 61 Gällstad

Tel. 0321-135 50

Fax. 0321-135 52

info@sableline.se

www.sableline.se

SABLELINE
Healthcare

SABLELINE
— S W E D E N —